


CONSIDERATIONS FOR INTERPRETING ANTIGEN TEST RESULTS IN NURSING HOMES



This algorithm should be used as a guide, but clinical decisions may deviate from this guide if indicated. Contextual factors including community incidence, characteristics of different antigen testing platforms, as well as availability and turnaround times of RT-PCR, further inform interpretation of antigen test results.

RT-PCR: reverse-transcriptase polymerase chain reaction

POC: point-of-care

HCP: healthcare personnel

Index case: a newly identified case of SARS-CoV-2 infection in a resident or HCP in a nursing home facility with no known infections of SARS-CoV-2 infection in the previous 14-day period.

COVID-19 outbreak response in a nursing home is triggered when one nursing home-onset SARS-CoV-2 infection in a resident or one HCP SARS-CoV-2 infection.

* Asymptomatic individuals who have recovered from SARS-CoV-2 infection in the past 3 months and live or work in a nursing home performing facility-wide testing do not need to be retested. If an individual has recovered from SARS-CoV-2 infection in the past 3 months and develops new symptoms suggestive of COVID-19, alternative diagnoses should be considered prior to retesting for SARS-CoV-2.

** Some antigen platforms have higher sensitivity when testing individuals within 5 days of symptom onset. Clinical discretion should be utilized to determine if retesting by RT-PCR is warranted.

[CMS recommendations](#) for testing asymptomatic HCP in facilities without a case
[CDC guidance on testing residents of nursing homes](#). [CDC guidance on testing HCP](#)

^ In discussion with the local health department, community incidence and time between antigen test and RT-PCR test can be utilized to interpret discordant results and determine when HCP can return to work.

^^ If an antigen test is presumptive negative in a facility with an outbreak, residents should be placed in transmission-based precautions or HCP should be allowed to continue working while monitoring for symptoms.

